

1D^2 genomic DNA sequencing for the MinION™ device using SQK-LSK308 (1/3)

Flow Cell Number
DNA Samples

Before start checklist <input type="checkbox"/> 1D^2 Sequencing Kit (SQK-LSK308) <input type="checkbox"/> Library Loading Beads Kit (EXP-LLB001) <input type="checkbox"/> SpotON Flow Cell FLO-MIN107 <input type="checkbox"/> Vortexer and Hula mixer <input type="checkbox"/> Pipettes P1000, P200, P100, P20, P10 and P2 <input type="checkbox"/> Pipette tips P2, P20, P100/200, P1000	<input type="checkbox"/> Approx 20 DNA LoBind Eppendorf 1.5 ml tubes <input type="checkbox"/> Freshly-prepared 70% EtOH <input type="checkbox"/> Covaris g-TUBE (optional) <input type="checkbox"/> AMPure beads resuspended and at RT <input type="checkbox"/> Magnet for bead separation <input type="checkbox"/> Nuclease-free water (NFW) <input type="checkbox"/> Laptop sleeptimer / update off <input type="checkbox"/> Microfuge <input type="checkbox"/> NEBNext FFPE Repair Mix (M66360) (optional)	<input type="checkbox"/> Run platform QC prior to library prep <input type="checkbox"/> NEB Blunt / TA Ligase Master Mix <input type="checkbox"/> 0.2 ml thin-walled PCR tubes <input type="checkbox"/> Heat block at 37 °C capable of taking 1.5 ml tubes <input type="checkbox"/> Thermal cycler at 20 °C and 65 °C <input type="checkbox"/> NEBNext Ultra II End-repair / dA-tailing Module <input type="checkbox"/> Update checks for MinKNOW™ completed <input type="checkbox"/> Check that the MinKNOW output directory is on the SSD
---	---	--

MASSFLOW	INSTRUCTIONS	NOTES / OBSERVATIONS	TIME / DATE
<p>DNA LoBind 45 µl</p> <p>Covaris g-TUBE</p> <p>DNA LoBind 45 µl</p> <p>17 µl</p> <p>DNA LoBind 62 µl</p> <p>62 µl</p> <p>DNA LoBind 124 µl</p> <p>Wash (2x) 200 µl</p> <p>46 µl</p> <p>DNA LoBind</p> <p>1 µl</p> <p>DNA LoBind 45 µl</p> <p>15 µl</p> <p>DNA LoBind 60 µl</p> <p>60 µl</p> <p>DNA LoBind 120 µl</p> <p>Wash (2x) 200 µl</p> <p>25 µl</p> <p>DNA LoBind</p> <p>1 µl</p> <p>DNA LoBind</p> <p>50 µl</p> <p>DNA LoBind 50 µl</p>	Genomic DNA 1-1.5 µg Genomic DNA in 45 µl Shear in a Covaris g-TUBE Eppendorf 5424; 6000 rpm for 8 kb fragments. 2 x 1 minute (invert tube to collect) Keep processing time under 15 minutes	This step is optional and can be omitted if long fragments are required.	
	DNA repair (optional) Add from NEBNext FFPE RepairMix <input type="checkbox"/> 8.5 µl NFW <input type="checkbox"/> 6.5 µl FFPE Repair Buffer <input type="checkbox"/> 2 µl FFPE Repair Mix		
	Mix by pipetting + spin down Incubate for 15 mins at 20 °C		
	<input type="checkbox"/> Add 62 µl resuspended AMPureXP beads to the end-prep reaction at room temperature and mix gently by pipetting <input type="checkbox"/> Incubate on rotator for 5 mins, spin down and pellet on magnet		
	<input type="checkbox"/> Keep on magnet wash 2x with 200 µl fresh 70% EtOH, do not disturb pellet <input type="checkbox"/> Spin down, replace on magnet, pipette of residual wash Briefly allow to dry		
	<input type="checkbox"/> Resuspend pellet in 46 µl NFW, incubate for 2 mins at RT <input type="checkbox"/> Pellet beads on magnet, remove eluate and transfer to fresh DNA LoBind tube		
	1 µl Qubit fluorometer - recovery aim < 1000 ng of material		
	Add from NEBNext Ultra II End-Repair / dA-tailing Module <input type="checkbox"/> 7 µl Ultra II End-Prep buffer <input type="checkbox"/> 3 µl Ultra II End-Prep enzyme mix <input type="checkbox"/> 5 µl NFW Mix by inversion + spin down Incubate for at 20 °C for 5 minutes and 65 °C for 5 minutes		
	<input type="checkbox"/> Add 60 µl resuspended AMPure XP beads at RT <input type="checkbox"/> Incubate on rotator for 5 minutes, spin down and pellet on magnet. Discard the supernatant. <input type="checkbox"/> Keep on magnet, wash 2x with 200 µl fresh 70% EtOH, do not disturb pellet <input type="checkbox"/> Briefly spin down, replace on magnet, pipette off residual wash. Briefly allow to dry <input type="checkbox"/> Resuspend pellet in 25 µl Nuclease-free water, incubate at RT for 2 minutes <input type="checkbox"/> Pellet beads on a magnet, remove eluate of End-Prepped DNA and transfer to fresh DNA LoBind		
	1 µl Qubit fluorometer – recovery aim about 700 ng of material		
	<input type="checkbox"/> Thaw the 1D^2 Adapter <input type="checkbox"/> Add 22.5 µl end prepped DNA (500 ng) <input type="checkbox"/> Add 2.5 µl 1D^2 Adapter <input type="checkbox"/> Add 25 µl Blunt/TA Ligase Master Mix Mix gently by inversion + spin down Incubate at RT for 10 minutes		

1D² genomic DNA sequencing for the MinION™ device using SQK-LSK308 (2/3)




Flow Cell Number

DNA Samples

MASSFLOW	INSTRUCTIONS	NOTES / OBSERVATIONS	TIME / DATE
	<p>Add 20 µl resuspended AMPure XP beads to the DNA sample</p> <ul style="list-style-type: none"> <input type="checkbox"/> Incubate on rotator for 10 minutes, spin down and pellet on magnet. Discard the supernatant. <input type="checkbox"/> Keep on magnet, wash 2x with 1 ml fresh 70% EtOH, do not disturb pellet <input type="checkbox"/> Briefly spin down, replace on magnet, pipette off residual wash. Briefly allow to dry <input type="checkbox"/> Resuspend pellet in 46 µl nuclease-free water, incubate at RT for 2 minutes <input type="checkbox"/> Pellet beads on a magnet, remove eluate of barcoded DNA and transfer to fresh DNA LoBind 		
	1 µl Qubit fluorometer		
	<p>Blunt/TA Ligase Master Mix</p> <p>When preparing the ligation reaction, mix by inversion between each addition</p> <ul style="list-style-type: none"> <input type="checkbox"/> 45 µl of 1D² adapted sample 5 µl <input type="checkbox"/> BAM <input type="checkbox"/> 50 µl Blunt/TA Ligase Master Mix <p>Mix gently by inversion and spin down</p> <p>Incubate at RT for 10 minutes</p>		
	<ul style="list-style-type: none"> <input type="checkbox"/> Vortex the AMPure XP beads to resuspend <input type="checkbox"/> Transfer 40 µl of beads to the adapter ligation reaction <input type="checkbox"/> Mix by pipetting, incubate at RT on a rotator mixer for 10 mins <input type="checkbox"/> Pellet beads on magnet and remove supernatant <input type="checkbox"/> Add 140 µl ABB to beads. Close tube lid, resuspend beads by flicking. Pellet beads on magnet and remove supernatant. Repeat. 		
	<p>Elute the adapted library (Pre-sequencing Mix)</p> <ul style="list-style-type: none"> <input type="checkbox"/> Resuspend pelleted beads in 15 µl ELB and incubate at RT for 10 minutes <input type="checkbox"/> Pellet beads on the magnet, remove the eluate and transfer to new DNA LoBind tube 		
	1 µl Qubit fluorometer – recovery aim about 200 ng of material		
<p>Before start checklist</p> <div> <input type="checkbox"/> MinION™ connected to computer with SpotON Flow Cell <input type="checkbox"/> Computer setup to run MinKNOW <input type="checkbox"/> Pre-sequencing Mix (PSM) at > 4 ng/µl </div> <div> <input type="checkbox"/> Run platform QC in parallel to library prep <input type="checkbox"/> Albacore downloaded <input type="checkbox"/> RBF and PSM on ice </div> <div> <input type="checkbox"/> Run Name set <input type="checkbox"/> NFW at RT <input type="checkbox"/> Platform QC completed </div>			
<p>Priming and loading the library</p>	<p>Prepare the MinION for sequencing protocol</p> <p>The platform QC should be run prior to library preparation beginning</p> <ul style="list-style-type: none"> <input type="checkbox"/> Assemble the MinION and MinION Flow Cell <input type="checkbox"/> Setup MinKNOW to run the Platform QC – name the run and start the protocol script – NC_Platform_QC.py <input type="checkbox"/> Allow the script to run to completion and the number of active pores are reported 		
	<p>Prime the flow cell ready for the library to be loaded when library preparation is complete</p> <p>Prepare priming buffer</p> <ul style="list-style-type: none"> <input type="checkbox"/> 480 µl RBF <input type="checkbox"/> 520 µl Nuclease-free water 		

1D² genomic DNA sequencing for the MinION™ device using SQK-LSK308 (3/3)

Flow Cell Number
DNA Samples

MASSFLOW	INSTRUCTIONS	NOTES / OBSERVATIONS	TIME / DATE
 <p>800 µl</p> <p>Priming Port</p> <p>5 minutes</p> <p>200 µl</p>	<p>Prime the flow cell using the priming port</p> <ul style="list-style-type: none"> <input type="checkbox"/> Open the priming port. Draw back a few µls of buffer to make sure there is continuous buffer flow from the priming port across the sensor array. <input type="checkbox"/> Load 800µl of the priming buffer. Wait 5 minutes <input type="checkbox"/> Gently lift the SpotON cover to make the SpotON port accessible <input type="checkbox"/> Load 200 µl of the priming buffer through the sample port 	<p>The aim of the draw back step is to remove any small air bubbles that have developed in the channel between the priming port and the sensor array window. This should be achieved by removing just a few µls of buffer before the air bubble is taken up into the pipette tip. Avoid removing large volumes that result in the sensor array no longer being covered by buffer.</p>	
<p>DNA LoBind</p> <p>35.0 µl</p> <p>25.5 µl</p> <p>12 µl</p> <p>2.5 µl</p>	<p>Prepare the library for loading</p> <ul style="list-style-type: none"> <input type="checkbox"/> 35.0 µl RBF kept at RT <input type="checkbox"/> 25.5 µl LLB kept at RT <input type="checkbox"/> 12.0 µl Adapted and tethered library <input type="checkbox"/> 2.5 µl Nuclease-free water <p>Mix gently by pipetting</p>	<p>RBF</p> <p>LLB</p>	
<p>DNA LoBind 75 µl</p> 	<p>Loading the prepared library</p> <ul style="list-style-type: none"> <input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON port in a dropwise fashion. <input type="checkbox"/> Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the sample port cover, making sure the bung enters the SpotON port <input type="checkbox"/> Close the sample port cover and replace the MinION lid. 		
	<p>Starting the sequencing script in MinKNOW</p> <ul style="list-style-type: none"> <input type="checkbox"/> Return to MinKNOW, name the run, select the NC_48Hr_Sequencing_Run_FLO-MIN107_SQK-LSK308.py and start using the start in the MinKNOW dialogue box. <input type="checkbox"/> MinKNOW will report the number of pores available for sequencing before data collection begins. These may differ from those reported in the Platform QC. <input type="checkbox"/> Allow the protocol to proceed until MinKNOW reports Finished Successfully. Use the Stop in the Control Panel to finish the protocol. <input type="checkbox"/> Open Albacore, and run the 1D² basecalling workflow 		