1D Genomic DNA sequencing for the MinION™ device using SQK-LSK108

Covaris g-TUBE (optional), NEBNext FFPE Repair Mix



Flow Cell Number: DNA Samples: Before start checklist Ligation 1D Sequencing Kit SQK-LSK108 ■ NEB Blunt/TA Ligase Master Mix Notebook sleeptimer/update off (AMX 1D in freezer until needed) ☐ AMPure beads Heat block for 1.5ml tubes at 37°C User adapt loading for (FLO-MIN104) Microfuge ☐ Thermal cycler at 20°C and 65°C NEBNext Ultra II End-repair/dA-tailing module Pre-chilled freezer pack for 0.2ml tubes at Vortexer and Hula mixer Magnetic rack for bead separation Pipettes P2, P20, P100/200, P1000 Nuclease-free water (NFW) (not supplied) ☐ DNA LoBind Eppendorf tubes 1.5ml Pipette tips P2, P20, P100/200, P1000 10mM Tris-HCl pH8.5 (if used) 0.2ml thin-walled PCR tubes 0.2ml thin-walled PCR tubes Latest versions Freshly prepared 70% EtOH Platform QC completed on flow cell of MinKNOW and Desktop Agent (check for updates) MinION SpotON FLO-MIN106 or FLO-MIN107 Additions to checklist for this protocol:

Covaris g-10BE (optional), NEBNext FFPE Repair Mix				
MASSFLOW	INSTRUCTIONS	NOTES/ OBSERVATIONS	TIME/DATE	
DNA LoBind 46 µl 1 µl Covaris g-TUBE 45 µl	Genomic DNA 1-1.5 µg Genomic DNA in 46 µl Shear in a Covaris g-TUBE Eppendorf 5424; 6000rpm for 8 kb fragments. 2 x 1 minute (invert tube to collect) Keep processing time under 15 minutes Analyse 1 µl for fragment size, quantity and quality (Agilent Bioanalyser)	This step is optional and can be omitted if long fragments are required.		
DNA LoBind 45 µl	DNA repair (optional) Add from NEBNext FFPE RepairMix (NEB M6630) 8.5 µl NFW 6.5 µl FFPE Repair Buffer 2 µl FFPE Repair Mix			
DNA LoBind 62 μl 62 μl	Mix by inversion + spin down Incubate for 15 mins at 20 °C			
DNA LoBind 124 µl	 □ Add 62 µl resuspended AMPureXP beads to end-prep reaction at RT □ Incubate on rotator for 5 mins, spin down and pellet on magnet Discard the supernatant at RT 			
Wash (2x) 200 µl	Keep on magnet wash 2x with 200 µl fresh 70% EtOH, do not disturb pellet Spin down, replace on magnet, pipette of residual wash Briefly allow to dry			
DNA LoBind	 ☐ Resuspend pellet in 46 µl NFW, incubate for 2 mins at RT Pellet ☐ beads on magnet, remove eluate and transfer to fresh DNA LoBind tube 			
1 µl	1 µl QuBit fluorimeter-recovery ≥ 1000 ng of material			
DNA LoBind 45 µl 15 µl DNA LoBind 60 µl	Add from NEBNext Ultra II End-Repair / dA-tailing Module 7 µI Ultra II End-Prep buffer 3 µI Ultra II End-Prep enzyme mix 5 µI DNA CS 3.6kb (DNA CS – positive control) Mix by inversion + spin down Incubate for at 20 °C for 5 minutes and 65 °C for 5 minutes	DCS		

nanoporetech.com Page 1/3

1D Genomic DNA sequencing for the MinION[™] device using SQK-LSK108



Flow Cell Number: DNA Samples:

MASSFLOW	INSTRUCTIONS	NOTES/ OBSERVATIONS	TIME/DATE	
DNA LoBind 120 μl Wash 2x 200 μl 31 μl	 Add 60 µl resuspended AMPure XP beads at RT Incubate on rotator for 5 minutes, spin down and pellet on magnet. Discard the supernatant. Keep on magnet, wash 2x with 200 µl fresh 70% EtOH, do not disturb pellet Briefly spin down, replace on magnet, pipette off residual wash. Briefly allow to dry Resuspend pellet in 31 µl Nuclease-free water, incubate at RT for 2 minutes Pellet beads on a magnet, remove eluate of End-Prepped DNA and transfer to fresh DNA LoBind 			
1 µl	1 µl QuBit fluorimeter - recovery aim about 700 ng of material			
DNA LoBind 30 µl 100 µl DNA LoBind 100 µl	NEB Blunt / TA Ligase Master Mix. When preparing the ligation reaction, mix by inversion between each addition. Check that the Master Mix is clear of any precipitate 30 µl end-prepped DNA 20 µl Adapter Mix 50 µl NEB Blunt / TA Master Mix Mix gently by inversion and spin down Incubate at RT for 10 minutes	AMX		
40 μl Wash 2x 140 μl	 □ Vortex AMPure XP beads to resuspend □ Transfer 40 µl of beads into the adapter ligation reaction □ Mix by pipetting, incubate at RT on a rotator mixer for 5 mins □ Pellet beads on magnet and remove supernatant □ Add 140µl ABB to beads. Close tube lid, resuspend beads by flicking. Pellet beads on magnet and remove supernatant. Repeat 	ABB		
15 µl	Elute the adapted library (Pre-sequencing Mix) Resuspend pelleted beads in 15 µl ELB and incubate at RT for 10 minutes Pellet beads on the magnet, remove the eluate and transfer to new DNA LoBind tube			
DNA LoBind 14 µl Store on ice	1 μl QuBit fluorimeter			
Before start checklist				
MinION™ connected to compu SpotOn Flow Cell	ter with Desktop Agent setup	RBF and library on ice		
Platform QC completed (car in parallel to library prep)	be done Run Name set	NFW at RT		
□ Computer setup to run MinKNOW □ Prepared library > 4 ng/μl				
Priming and loading the library MinKNOW	Prepare the MinION for sequencing protocol The platform QC should be run prior to library preparation beginning Assemble the MinION and MinION Flow Cell Setup MinKNOW to run the Platform QC – name the run and start the protocol script – NC_Platform_QC.py Allow the script to run to completion and the number of active pores are reported			

nanoporetech.com Page 2/3

1D Genomic DNA sequencing for the MinION™ device using SQK-LSK108



Flow Cell Number: DNA Samples:

MASSFLOW	INSTRUCTIONS	NOTES/ OBSERVATIONS TIME/DATE
Sample cover Activator	Prime the flow cell ready for the library to be loaded when library preparation is complete Prepare priming buffer 480 µl RBF 520 µl Nuclease-free water	
200 µl	Prime the flow cell Open the sample port. Draw back a few µls of buffer to make sure there is continuous buffer flow from the sample port across the sensor array. Load 800µl of the priming buffer. Wait 5 minutes Gently lift the activator to make the SpotON port accessible Load 200µl of the priming buffer through the sample port	RBF
DNA LoBind 200 µl 35.0 µl 25.5 µl 2.5 µl	Prepare the library for loading 35.0 µl RBF kept at RT 25.5 µl LLB kept at RT 12.0 µl Adapted and tethered library 2.5 µl Nuclease-free water Mix gently by pipetting	RBF (LLB)
DNA LoBind 75 µl	Loading the prepared library Add 75µl of sample to the flow cell via the SpotON port in a dropwise fashion. Ensure each drop flows into the port beforeadding the next. Gently replace the activator, making sure the bung enters the SpotON port Close the sample port cover and replace the MinION li	
MinKNOW	Starting the sequencing script in MinKNOW and the workflow in the Desktop Agent Return to MinKNOW, name the run, select the NC_48Hr_Sequencing_Run_FLO_MIN106_SQK-LSK108 _plus_Basecaller.py for live basecalling) using the start in the MinKNOW dialogue box MinKNOW will report the number of pores available for sequencing before data collection begins. These may differ from those reported in the Platform QC. Allow the protocol to proceed until MinKNOW reports Finished Successfully. Use the Stop in the Control Panel to finish the protocol. Quit the Desktop Agent, close down MinKNOW and disconnect the MinION After sequencing	
	If using Albacore for local basecalling please refer to the instructions in Albacore basecalling software	
After sequencing check	list	
Store washed flow cell at 4°C or the returns form in the Nanopore		Navigate to www.metrichor.com to review the full sequencing report
Store MinION at RT		

nanoporetech.com Page 3/3