

# MinION workshop - Getting to know your data

```
pandoc brachy.md --smart --standalone --bibliography test.bib -o brachy.pdf
```

## Section Goals

- Understand the initial data that is coming from your MinION run
- What is contained within the **fast5** datafiles
- What we can, and cannot, do with the data at this point

## General pipeline so far:

- The MinKNOW software on the sequencing laptop facilitates the collection and parsing of data from the MinION flowcell.
- This ‘raw’ data is voltage traces that are measured at 5khz (5 thousands times per second).
- The voltage across the membrane changes as molecules enter, move through, and exit each sequencing pore.
- This entire measurement is called an ‘event’ which will be the collected signal measurements for that molecule.
- The software recognizes this initial voltage change and begins writing out the trace to a new **fast5** file.
- There is **one file per DNA molecule** that passes through each pore. Therefore many files are created!
- Hence, the requirement for high quality transfer speeds on your laptop to facilitate this (SSDs)
- This is a fundamental shift from Illumina platforms in which most analysis is based off of fewer, but much larger, files.

## The initial data is NOT basecalled

- At this point, we do not have nucleotide information, just voltage traces.
- We can look at this information within the fast5 files that are created.

## What are **fast5** files?

- Each file contains signal trace for a single DNA molecule
- **fast5** files are stored in the **HDF5** file format.
- **HDF5** acts as a larger organizational framework (box) in which different types of data can be stored in a hierarchical way.
- Like directories, sub-directories, and files/data on your computer.
- The **HDF5** format is quite flexible in being able to contain a wide variety of data types internally (strings,int,float,arrays, etc) within the single ‘**HDF5**’ format.
- The **fast5** files are therefore able to not only hold the signal trace, but also a wide variety of metadata related to how it was created, timestamps, software versions, etc.
- The format allows for new items to be added to it over time, creating a single container from raw information onwards
- We will see this shortly.
- We can investigate any **fast5** file using an **HDF5** ‘viewer’ which will show the internal structure and datasets in an easy-to-examine GUI.
- The group that created the **HDF5** format [provides a java-based viewing software to download](#):
- [HDFView for windows](#)
- [HDFView for OSX](#)
- [HDFView for Linux](#)
- There are also [command-line tools to viewing HDF5 files for those interested](#)

## Where are the **fast5** files located?

- As sequencing begins, data beings to be filled in to the **./reads** directory within your MinKNOW install folder:

As of MinKNOW v1.6.11:

Windows: OSX: **/Library/MinKNOW/data/reads/** Linux:

## Let’s take a look at a **fast5** file




Figure 1: screenshot from hdfview